



**MELAMINE IN MILK PRODUCTS FOR HPLC ANALYSIS USING:  
60mg 3mL HyperSep Retain-CX Extraction column  
(Part Number: 60107-303)**

**1. PREPARE SAMPLE**

Weigh 5g of milk powder (or measure 10mL of milk) into a 250mL flask.  
Add 50mL of 1% Trichloroacetic acid (TCA).  
Mix / vortex.

Add 2mL of 2% lead acetate / water solution into the mixture then sonicate for 20 minutes.  
Transfer part of the final mixture into a 10mL centrifuge tube. Centrifuge for 10 minutes at 8,000 rpm.

**2. CONDITION RETAIN-CX EXTRACTION COLUMN**

3mL of methanol  
3mL of DI water

**3. APPLY SAMPLE**

Load 6mL of the sample extract onto the Retain-CX SPE product at 1-2mL/minute.

**4. WASH COLUMN**

3mL of DI water  
3mL of methanol  
Dry column (5 minutes at > 10 inches Hg).

**5. ELUTE MELAMINE**

5mL of 5% ammonia / methanol  
Collect eluate at 1 to 2mL/min.

**6. DRY ELUATE**

Evaporate to dryness at < 50°C using nitrogen.  
Reconstitute sample using 1mL of mobile phase