

# Determination of Melamine in Powdered Milk by LC-MS/MS Using a Core Enhanced Technology Solid Core HPLC Column

Kimberly Phipps, Thermo Fisher Scientific, Runcorn, Cheshire, UK

## Key Words

Melamine, milk, HyperSep, Retain-CX, Accucore HILIC


## Abstract

In recent years melamine has been found in adulterated milk and milk based products, especially in infant formula milk, causing thousands of children to become ill. Limits have now been set for the amount of melamine that can be present in foodstuffs; for infant formula milk this level is zero and for other food stuffs there is a 2.5 ppm maximum. An LC/MS method for melamine in infant formula milk has been developed. Using Thermo Scientific™ HyperSep™ Retain-CX solid phase extraction cartridges sample preparation is fast and efficient giving an average recovery of 121%. The Thermo Scientific Accucore™ HILIC column provided a fast run time of 2 minutes. The dynamic range was linear between 10 and 1000 ng/g with a correlation coefficient ( $r^2$ ) of 0.9913 and accuracies of +/- 20% for standards.

## Introduction

Melamine is a compound commonly used as a fire retardant in plastics or as a fertilizer due to its nitrogen rich structure. As melamine is nitrogen rich it has been used to fraudulently increase the detected protein level in milk so that the milk can command a higher price. The ingestion of melamine can cause renal complications including kidney and bladder stones and even bladder cancer. In 2008, (13000) children were taken ill in China as a result of drinking infant formula milk containing melamine. Limits have now been set for the amount of melamine that can be present in foodstuffs, for infant formula milk this level is zero (although a level below 1 ng/g does not present a concern) and for other food stuffs is 2.5 ppm. It is therefore necessary to have a fast and sensitive method that can be conducted quickly and gives reproducible results.

The method described in this application note allows the testing of infant formula milk for the presence of melamine. Using HyperSep Retain-CX for the extraction of melamine followed by analysis using an Accucore HILIC column, a fast two minute run time was achieved.



## Experimental Details

| Consumables  | Part Number |
|--|-------------|
| Fisher Scientific LCMS grade methanol                        | M/4062/17   |
| Fisher Scientific LCMS grade water                           | W/0112/17   |
| Fisher Scientific LCMS grade acetonitrile                    | 51101       |
| Fisher Scientific analytical reagent grade formic acid 100 % | F/1900/PB08 |
| Fisher Scientific HPLC grade ammonia solution                | A/3295/PB05 |
| Fisher Scientific analytical reagent grade acetic acid       | A/0415/07   |
| Melamine, Sigma Aldrich                                      | 30130TR     |
| Infant Formula Milk – Stage 2                                |             |

| <b>Sample Handling Equipment</b>                |  | <b>Part Number</b> |
|---|--|--------------------|
| Liquid handling hardware:                       |  |                    |
| Thermo Scientific FinnPipette (100-1000 µL)     |  | 642090             |
| Thermo Scientific FinnPipette (10-100 µL)       |  | 4642070            |
| Thermo Scientific FinnPipette (1-10 µL)         |  | 4642040            |
| <b>SPE Hardware</b>                             |  | <b>Part Number</b> |
| Thermo Scientific Ultra Vap                     |  | CLS-229070         |
| Thermo Scientific HyperSep Glass Block Manifold |  | 60104-232          |
| <b>SPE Cartridge</b>                            |  | <b>Part Number</b> |
| HyperSep Retain-CX, 60 mg/3 mL                  |  | 60107-303          |
| <b>Vials and Closures</b>                       |  | <b>Part Number</b> |
| Thermo Scientific Chromacol storage vial        |  | 22-SV-CP           |
| Thermo Scientific premium vial                  |  | 60180-600          |
| <b>Sample Preparation</b>                       |  |                    |
| Compound:                                       | Melamine   |                    |
| Matrix:   | Milk powder  |                    |
| Weigh   | 1 g of milk powder. Dissolve in 4 mL deionised water. Mix well   |                    |
|   | Spike in appropriate amount of melamine calibrator.<br>Add 6 mL of 2.5% formic acid (aq). Mix well.<br>Sonicate for 10 minutes and centrifuge for 20 minutes at 2000 rpm         |                    |
| SPE cartridge type:                             | HyperSep Retain-CX 60 mg/3 mL  |                    |
| Conditioning stage:                             | Load 3 mL of methanol onto the SPE cartridge under vacuum followed by 3 mL of water  |                    |
| Application stage:                              | Load all supernatant onto the SPE cartridge, trying not to disturb the top layer   |                    |
| Washing stage:                                  | Load 3 mL of water followed by 3 mL of methanol onto the SPE cartridge   |                    |
| Elution stage:                                  | Load 2 x 1 mL 5% ammonia in methanol onto the SPE cartridge  |                    |
| Additional stage:                               | Evaporate the collected sample to dryness under nitrogen.<br>Reconstitute in 1 mL 90:10 acetonitrile: 50 mM ammonium acetate pH 5. Gently vortex mix and sonicate for 5 minutes. |                    |
| <b>Separation Conditions</b>                    |  | <b>Part Number</b> |
| Instrumentation:                                | Thermo Scientific Accela™ 600 pump,<br>Thermo Scientific CTC autosampler   |                    |
| Column:   | Accucore HILIC 2.6 µm, 100 x 2.1 mm  | 17526-102130       |
| Mobile phase A:                                 | 50 mM ammonium acetate adjusted to pH 5 with acetic acid   |                    |
| Mobile phase B:                                 | Acetonitrile   |                    |
| Isocratic:                                      | (10:90)  |                    |
| Flow rate:                                      | 1 mL/min   |                    |
| Run time:                                       | 2 minutes  |                    |
| Column temperature:                             | 40 °C  |                    |
| Injection details:                              | Inject 10 µL   |                    |
| Loop Size:                                      | 100 µL   |                    |
| Injection wash solvent:                         | Acetonitrile: water (50:50)  |                    |
| Preparation of mobile phase A:                  | Weigh 3.854 g of ammonium acetate and dilute into 1 L of water.<br>Adjust to pH 5 with acetic acid   |                    |

**MS/MS Conditions**

|                         |                               |
|-------------------------|-------------------------------|
| Instrumentation:        | Thermo Scientific TSQ Vantage |
| MS/MS run time:         | 2 minutes                     |
| Ionization conditions:  | APCI Positive                 |
| Discharge current:      | 4 eV                          |
| Vaporizer temperature:  | 350 °C                        |
| Sheath gas pressure:    | 50 arbitrary units            |
| Auxiliary gas pressure: | 15 arbitrary units            |
| Capillary temperature:  | 300 °C                        |
| Collision pressure:     | 1.2 mTorr                     |

**Ions Monitored**

| Precursor (m/z) | Product (m/z) | Collision Energy (eV) | S-lens (Arb) |
|-----------------|---------------|-----------------------|--------------|
| 127.1           | 68.213        | 32                    | 32           |
| 127.1           | 85.171        | 17                    | 44           |

|                |                   |
|----------------|-------------------|
| Scan time:     | 0.1 s             |
| Q1 peak width: | 0.20 m/z          |
| Q3 peak width: | 0.70 m/z          |
| Chrom filter:  | 5 arbitrary units |

**Solutions**

|                 |   |
|-----------------|---|
| Melamine stock: | Weigh 1 mg and dissolve in 1 mL of 50 mM ammonium acetate pH 5. |
|-----------------|---|

**Data Processing**

|                        |                              |
|------------------------|------------------------------|
| Software:              | Thermo Scientific LCQuan 2.6 |
| Integration algorithm: | ICIS                         |

| Parameter                           | Value                                |
|-------------------------------------|--------------------------------------|
| MRM                                 | 127.1 (68.208-68.218, 85.166-85.176) |
| Retention Time (minutes)            | 0.59                                 |
| Retention Window (seconds)          | 15                                   |
| Smoothing (arbitrary units)         | 1                                    |
| Baseline Window (counts)            | 75                                   |
| Area Noise Factor (arbitrary units) | 28                                   |
| Peak Noise Factor (arbitrary units) | 10                                   |
| Constrain Peak Width                | No                                   |
| Component Type                      | Target Compound                      |
| Calibration Curve Type              | Linear                               |
| Response                            | Area                                 |
| Origin                              | Ignore                               |
| Weighting                           | 1/X <sup>2</sup>                     |

## Results

### Chromatography

The Accucore HILIC column gave excellent peak shape. The chromatography of Melamine can be seen in figure 1 and the chromatography of a milk blank can be seen in figure 2.




Figure 1: Selected Reaction Monitoring chromatogram of melamine at 500 ng/mL




Figure 2: Selected Reaction Monitoring chromatogram of milk blank

### Linearity

Standards of melamine extracted from spiked milk powder gave a linear calibration curve over the dynamic range of 10 to 1000 ng/mL with an  $r^2$  of 0.9928 (figure 3 and table 1).




Figure 3: Linearity over the dynamic range of 10 ng/mL to 1000 ng/mL

| Sample       | Melamine Response (peak area) | Specified Concentration (ng/g) | Calculated Concentration (ng/g) | % Diff |
|--------------|-------------------------------|--------------------------------|---------------------------------|--------|
| Extracted S1 | 7770                          | 10                             | 10                              | 1      |
| Extracted S2 | 14549                         | 25                             | 24                              | -4     |
| Extracted S3 | 28124                         | 50                             | 52                              | 4      |
| Extracted S4 | 50881                         | 100                            | 98                              | -1     |
| Extracted S5 | 129277                        | 250                            | 260                             | 4      |
| Extracted S6 | 205400                        | 500                            | 420                             | -17    |
| Extracted S7 | 411165                        | 750                            | 840                             | 12     |
| Extracted S8 | 492232                        | 1000                           | 1000                            | 1      |

Table 1: Accuracy of extracted standards

### Recovery

An average Recovery of 121% was observed.

### Conclusion

Retain CX SPE cartridges and Accucore HILIC HPLC columns allow for a simple extraction and rapid quantification melamine from milk powder.

Accucore HILIC gives a fast runtime of 2 min with excellent peak shape for melamine.

An LLOQ of 10 ng/mL was achieved with an extraction recovery of above 100%.

[thermoscientific.com/chromatography](http://thermoscientific.com/chromatography)

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**USA and Canada** +1 800 332 3331  
**France** +33 (0)1 60 92 48 34  
**Germany** +49 (0) 2423 9431 20 or 21  
**United Kingdom** +44 (0)1928 534110  
**Japan** +81 3 5826 1615

**China** +86 21 68654588 +86 10 84193588  
+86 20 83145199 800 810 5118  
**India** +91 22 6742 9494 +91 27 1766 2352  
**Australia** 1 300 735 292 (free call domestic)  
**New Zealand** 0800 933 966 (free call domestic)  
**All Other Enquiries** +44 (0) 1928 534 050

**Technical Support**  
**North America** +1 800 332 3331  
**Outside North America** +44 (0) 1928 534 440

**Thermo**  
**SCIENTIFIC**

Part of Thermo Fisher Scientific