


Innovative Solutions for High density sample handling

Thermo Scientific Well Plate Sealing Tapes for Chromatography

Do you face any of these issues with your tape?

- Your tape is not resistant against organic solvents.
- The tape is not resealing after puncture.
- The needle lifts the tape or plate.
- The adhesive contaminates your sample.



If so, we are able to offer the right solution to improve your flexibility and reproducibility working with well plates

Our 96-well plate compatible tapes for chromatography offer:

- Chromatography optimized material
- Pre-slitted but high resealability
- Polar organic solvent stability
- Acid stability
- No glue or adhesive in the cavity area

Description	Color	Material	Pre-slitted	Cat. No.	Pack of:
Sealing Tape with three layers (no adhesive at cavity area)	Clear	PET/Silicone/PET	Yes	60180-M146	100
Sealing Tape with three layers (no adhesive at cavity area)	Clear	PET/Silicone/PET	Yes	60180-M176	25

For more information about these products, visit
thermoscientific.com/webseal

©2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

USA and Canada +1 800 332 3331
France +33 (0)1 60 92 48 34
Germany +49 (0) 2423 9431 452
United Kingdom +44 (0)1928 534110
Japan +81 3 5826 1615

China +86 21 68654588 +86 10 84193588
+86 20 83145199 800 810 5118
India +91 22 6742 9494 +91 27 1766 2352
Australia 1 300 735 292 (free call domestic)
New Zealand 0800 933 966 (free call domestic)
All Other Enquiries +44 (0) 1928 534 050

Technical Support
North America +1 800 332 3331
Outside North America +44 (0) 1928 534 440

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand